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## Centrifuge technology revolutionises adhesion testing

Prof. Dr. Dietmar Lerche, Dipl. - Ing. Uwe Dietz - LUM

## Introduction

Composite materials, manufactured by bonding, layering, and coatings fulfil a great variety of functions in a wide field of applications. Bonding – as an emerging joining technology is applied in lightweight construction, in automotive and electronics industry, in optics, composite materials and many other fields. In addition, surface modification and coatings

provide special functions like corrosion resistance, special optical or electric properties or just a good appearance due to decorative coating. The strength of both, adhesive bonds and coatings, under the respective environmental conditions, is an essential for its functionality. It needs to be analysed both in design and quality control. However, bonding strength is commonly tested by tensile strength instru-

ments designed to handle only one sample. To achieve good reproducibility, special clamping aids and alignment constructions are necessary and joints of certain materials are difficult to test.

In case of coatings and multilayered materials common test procedures (e.g. peel test) allow only semi-quantitative or qualitative evaluation. This article describes a completely new quantitative test technology based on an-

alytical centrifugation. The new analytical instrument Adhesion Analyser LUMiFrac® uses the centrifugal force as test load and allows to determine bonding and adhesive strength of up to 8 samples simultaneously.

## State-of-the-art testing procedures

At present, adhesive or bonding strength is mostly deter-

mined either by pulling the two parts of the joint apart (DIN EN 15870, pull test, see fig. 1[1]), by scratching a coated surface according to a special lattice design (ISO 2409, cross-cut test [2]), or by peeling a coating off the substrate (ISO 8510 [3,4], Fig. 2). The tensile test is the only procedure that delivers the critical tensile strength, which causes the failure of joint or coating. But to achieve this, specially designed test

Fig. 2 Test setup di un peel a 90° secondo ISO 8510-1:1990 [3]  
 Test setup of a 90° peel test according to ISO 8510-1:1990 [3]

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*specimen must be used to clamp the two parts and care has to be taken to have a good alignment to avoid shear forces. This test method is only designed for single sample testing. Series of experiments to get a significant ranking of bonding strength between different material combinations, different adhesives, and/or curing conditions are therefore very time consuming.*

### Principle of the centrifuge technology

*In the centrifuge the radial directed centrifugal force provides the tensile force on a sample (Fig. 3).*

*The centrifugal force acting on the layer of adhesive is given by*  
 $FC = m \omega^2 r = m (2 \text{ RPM}/60)r$   
 (eq. 1)

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mass of the rotating test stamp,  $\omega$  the rotational angular velocity and  $r$  the radial distance of the centre of mass of the test stamp. With the effective cross-area of the glued parts,  $A$ , the tensile stress is given by eq. 2:

$$\sigma = F / A \quad (\text{eq. 2})$$

The loading force can easily be varied by the rotational speeds up to 13,000 RPM and the mass of the test stamp.

A typical test specimen consists of substrate1/adhesive/substrate2. Substrate1 itself acts as the test stamp and is glued by means of the adhesive to substrate2. The latter is supported by a cylindrical sleeve and

the only geometrical constraint consists in being in area larger as the cross section of the sleeve. The sleeve guides the stamp and this assembly is easily inserted into the failure detection modul (FDM). 8 of these modules, equipped with the same kind or different samples, are slid into the 8 detection positions of the rotor. This design assures a fixed position of each assembly without any clamping that may cause unwanted shear forces. Test stamps are available in different materials to cover a wide range of material combinations on the one hand side, and to enable a wide range of loads posed

upon the test specimen on the other hand side. With different test stamp masses, centrifugal forces from 0.1 N to 6.5 kN can be realised. In addition, strength can be varied by a factor of 4 by different cross areas of the test stamp (cp. eq.2). The Adhesion Analyser LUMiFrac comes with the approved database-based software SEPView®, which has already been used in other analytical centrifuges like LUMiSizer [7, 8]. It allows programming special SOPs (e.g. velocity of force increase, alternating loads, temperature), controls the measurement, performs the data acquisition and data

management as well as data analysis and documentation. A typical measurement starts with the input of the sample names, the selection of the type of test stamps from the library and measurement like RPM-time-progress. During the centrifugal test the rotational velocity is ramped up according to program settings. If the critical tensile strength of the test specimen is reached, the joint fails, coating fractures or delamination takes place and correspondingly the test stamp (cp. Fig. 3) moves slightly outward. This relocation is detected in situ by the corresponding FDM and a signal is communicated to

the software. Based on the continuously recorded RPM and the type of test stamp SEPView automatically calculates the critical force/strength of failure for each of the 8 samples according to eqs.1 and 2. Results are immediately displayed as a bar chart on the monitor (Fig. 4) and may be saved and printed together with other experimental data. LUMiFrac allows, on the one hand, the measurement of samples of the same type of materials combination or the comparison of samples regarding different substrates, adhesives, coatings and variations of processing. Despite the multisample approach,

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*the benchtop instrument is very compact (Fig. 5), and demands only a minimum hands-on-time. Software settings, loading the rotor, determination of the critical bonding strength/stress*

*according to DIN EN 15870 need less than 5 minutes. Additionally, reproducibility is not only likely to be improved as no influence of shear forces (clamping of pull test) is given, but also because all*

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8 samples are measured at the same time and the same environmental conditions.

## Application examples and results

The functionality of the measurement system has been verified with different joints, regarding material combination, adhesives and curing processes. Fig. 6 reveals the

influence of mixing ratio and curing temperature on the strength of 2-C-Epoxy resins bonding copper and construction steel S235.

It is easily to read off that higher adhesive strength can be achieved by modifying the mixing ratios as well as by curing at elevated temperatures. Parallel testing of bonds prepared from different material combinations and adhesives allows for direct comparison, as shown in fig. 7.

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Critical failure force (left) and corresponding strength (right) of differently coated glasses (mirrors) from the optical industry (see [10])

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*Professor Dr. sc. nat. D. Lerche is the founder and chairman of LUM (Labor, Umweltdiagnostik & Medizintechnik) GmbH based in Berlin, Germany. He attended Lomonossow University in Moscow, Russia, and received his PhD in biophysics from the well-known Math and Natural Science Faculty at Humboldt University, Berlin. In 1994 he founded the independent LUM GmbH with the goal to provide innovative and unique instrumentation to advance the fields of suspension and emulsion sciences.*

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*Uwe Rietz (Dipl.-Ing. (FH)) is an Development Engineer for LUM in Berlin, Germany. He obtained his degree in Automotive Engineering from the HTW Berlin (University of Applied Sciences).*

adhesive exhibits independently on the material combination the best tensile strength of more than 50 N/mm<sup>2</sup> (50 MPa). On the other hand, tensile strength for CY and PU is much lower and its ranking depends on the material combination.

Besides different materials, adhesives and different temperatures also different pre-treatment processes can be evaluated in one experiment. If the strength of a coating is to be tested, one has to make sure that the test stamp fixed by an adhesive to the surface of the coating reveals higher values of tensile strength than the coated sample or multi-layered material itself. In other words, the FDM

will detect any form of material failure that occurs for any critical force lower than that generated by the LUMiFrac (6.5 kN). This is not the case for common peel test (ISO 8510 [3,4]) or tape-cross-hatch test (ISO 9211-4 [9]), which are limited due to adhesive properties of the tapes.

Fig. 8 shows results obtained for a set of 9 different metal coated glasses. All of these samples passed the peel test according to ISO 9211-4. During the centrifugal test 5 different types of failure and fracture pattern occurred. Optical inspection revealed delamination fracture, adhesive and cohesive fracture, interlayer delamination, partial delamination fracture

and substrate fracture, listed in the order from the lowest to the highest stress that has been exerted on the samples to achieve these types of failure.

For details refer to [10].

### Conclusions

A new multisample measuring principle to determine bonding or adhesive strength of a wide range of bonded materials, coated materials and composite materials is described. In contrast to common pull and peel tests the Adhesion Analyser LUMiFrac® employs the radially directed centrifugal force to stress up to 8 different samples simultaneously. The LUMiFrac® is a

compact and easy to handle test apparatus, where test samples in contrast to e.g. common pull test need not to be clamped. The usually prepared test specimen is plugged into a special designed sleeve, which supports and guides the test stamp, slid into the failure detection module and placed into the rotor. Increasing or alternating test forces are applied by a programmed rotational velocity and the critical failure force is detected during centrifugation and communicated to the PC-software SEPView. Experimental data of each sample are automatically analysed, displayed online and saved. In comparison to a conventional tensile testing machine (e.g. DIN EN

15870, pull test), the total test time (software settings, putting in samples, test procedure and calculation of failure force/strength) for 8 samples is reduced to less than 85% (5 minutes). Experimental results regarding different adhesives, substrate combinations and metal coated glasses demonstrate the broad applicability of the new test procedure.

### References

- [1] DIN EN 15870: Adhesives – Determination of tensile strength of butt joints
- [2] ISO 2409:2007: Paints and varnishes – Cross-cut test (German version DIN EN ISO 2409:2007)
- [3] ISO 8510-1: 1990: Ad-

- hesives; peel test for a flexible-bonded-to-rigid test specimen assembly - part 1: 90° peel (German version DIN EN ISO 8510-1: 2006: Adhesives; peel test for a flexible-bonded-to-rigid test specimen assembly - part 2 180° peel (German version DIN EN ISO 8510-2:2010)
- [5] U.Rietz, S. Hielscher, U. Beck, D. Lerche; “A new centrifuge technology for efficient and reliable testing of composite strength of bonded joints and coated surfaces”; Congress-ACHEMA 2012; Session: Materials and Materials testing, June 20th, Frankfurt
- [6] U.Rietz, S. Hielscher, U. Beck, D. Lerche; Centrifuge technology revolutionises adhesion testing. Dispersion

- Letters, 4, 2013, 1-9, <http://www.dispersion-letters.com/>
- [7] D. Lerche, T. Sobisch; “Direct and Accelerated Characterization of Formulation Stability, J. of Dispersion Science and Technology, 32:12, 2011 1799-1811; <http://dx.doi.org/10.1080/01932691.2011.616365>
- [8] <http://www.lum-gmbh.com/dispersion-and-stability-analyser.html>
- [9] ISO 9211-4:2012: Optics and photonics – Optical coatings - Part 4: Specific test methods
- [10] U. Beck, G. Reiners, M. Weise, D. Lerche, U. Rietz, H. Niederwald; “Quantitative Adhesion Testing of PVD-Coatings by means of Centrifuge Technology”, Surface & Coatings Technology 205 (2011)182-186.